

 UNIVET



# LOUPES

Magnification Systems

# HEADGEAR

## SIZE

Adjustable Fitting

Magnification	 Galilean	 Prismatic
2.0X		
2.5X		
3.0X		
3.5X		•
4.5X *		•
6.0X *		•

\* Magnification available only with Headgear & Air-X systems

## ACCESSORIES INCLUDED

Magnetic headlight adapter

## FEATURES

- Flexible structure with no pressure points
- Adjustment of the vertical axis and declination angle
- Interchangeable sweatband sponges
- Can be worn over prescription eyewear

*Struttura flessibile, senza punti di pressione  
Regolazione di asse verticale e angolo di inclinazione  
Spugne parasudore intercambiabili  
Sovrapponibile agli occhiali da vista*

## FRAME

Details



Head  
gear

Available Prismatic

[www.univetloupes.com](http://www.univetloupes.com)

 UNIVET



UNIVET LOUPES Spa Via Giovanni Prati, 87 - 25086 Rezzato (BS) Italy +39 030 2499411 · [info@univetloupes.com](mailto:info@univetloupes.com)